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1. Product's name:

1.1. Technical Name: Plate for Making Molds

1.2. Commercial name: Memoflex PETG

1.3. Commercial Model: Memoflex PETG board 0.30 mm, 0.40 mm, 0.50 mm, 0.60 mm, 0.75 mm, 0.76 mm, 0.80 mm, 1.00 mm.

1.4. Anvisa Notification: 10331430036

2. Manufacturer:

Aditek do Brasil S.A.

Address: Cesário Mota, 14

CEP: 14140-000 - Cravinhos - SP - Brazil

Phone: +55 (16) 3951-9355 / 0800-0159520

CNPJ 64.602.097/0001-95

Technical Responsible: José Arantes Nogueira Neto - CREA: 601037885

www.aditek.com.br

3. Certifying entity:

Not applicable.

4. Usage Indication:

Memoflex plates are indicated for the manufacture of dental trays and according to their characteristics and thickness, for different dental applications such as orthodontic aligners, occlusion plate, tray for indirect bonding of brackets and tooth whitening plate.

5. Users:

The use of plates for making trays is restricted to professionals in the field of dentistry. The use by people outside this context may cause harm to the patient. The user undertakes to verify his responsibility for the suitability of the product for the intended use, especially in the case of fields of use not expressly indicated in the instructions for use.

6. Performance Features:

- Resistant to shocks and drops, even in reduced thicknesses;
- Excellent shape maintenance after vacuum thermoforming (dimensional stability);
- Inert to the action of most chemicals;
- Excellent transparency;
- It practically does not cause allergic reactions or tissue irritation on contact;
- Hygroscopic, that is, it easily absorbs moisture.

7. Clinical Benefits:

Thermoformed polymer: Fidelity in the reproduction of the arcades.

Excellent translucency: Aesthetics during use.

8. Warnings:

 Aditek products are designed and manufactured for single use only. Its reuse or reprocessing can cause cross-infection and decrease in efficiency.

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 If the packaging has been violated, do not use the product, avoiding risks related to the transmission of microorganisms or the possible falsification of products.

 Store the product in a place at room temperature, dry and away from sunlight, keeping it in its own packaging, so as not to cause changes in its shape.

This product must be transported in its own packaging, must not be folded or crumpled and protected from moisture.

The product is supplied in a non-sterile condition.

The professional must use good clinical practices for the disposal of medical products according to current health regulations.

Aditek do Brasil S.A. is not responsible for the results obtained by applying the product without proper precautions or failure to observe the warnings and instructions for use.

If serious incidents occur with the product, Aditek do Brasil and the competent authority at the national level must be reported.

9. Contraindications:

Not applicable.

10. Residual Risks

Not applicable.

11. Unwanted Side Effects

During the laminating process, do not neglect the laminating equipment. The material of the plates heats up and softens quickly and can become unusable and even cause damage to the plasticizing equipment if it softens too much.

When working with any model of plate for making trays, during manufacturing and finishing the tray, always wear a mask for your safety.

People who are not used to vacuum thermoforming procedures must read the instructions for use provided by the manufacturer before using the equipment.

The responsibility for the clinical treatment, technique for the construction of the tray, side effects and the expected action when using the trays generated by the Memoflex plates, are the sole responsibility of the dentist and/or the professional responsible for their use.

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12. Composition:

Polyethylene Glycol Terephthalate (PETG)
Chemical composition limits (%)
100%

13. Validity:

 2 years stored in its own packaging and closed.

14. Symbols and Meanings:

	Manufacturer		Product reference code
	Consult instructions for use		Expiration date
	Heads up		Do not reuse the product
	Do not use if package is damaged or tampered with		Lot number
	Keep away from sunlight		

15. Information to be given to the Customer:

The vacuum-sealed, metallized plastic packaging protects the boards from moisture and contaminants. The plates made of PETG are highly hygroscopic (high moisture absorption capacity). After opening the vacuum packaging, the plates that are not used immediately must remain stored in the same packaging so as not to absorb moisture and not present numerous small bubbles during heating in the laminator.

16. Usage mode:

All tray making plates are wrapped with insulating film and can be processed in a pressure or vacuum thermoforming unit.

If processed in a pressure thermoforming machine, this extremely thin film must be thermoformed together with the plate, on the model side, ensuring an improved transparency of the plate, without the need for additional insulation of the model surface. After thermoforming, the film must be removed.

If processed in a vacuum thermoformer, the insulating film must be removed prior to thermoforming.

Cleaning of the product should only be carried out after finishing the plate.

The heating and cooling time setting depends on the thickness of the thermoforming sheet and can be programmed directly on the thermoforming unit.

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Plasticizing Ideal Point:

Due to the existence of several types of plates with different thicknesses, the ideal point of plasticization is identified by changes in the plate, thus occurring different heating times. In general, the sweet spot is observed by changing the texture and surface color (gloss). For less translucent boards, the board will be shiny and transparent. For thicker or opaque slabs, the ideal point is to drain the slab by 10 to 12 mm.

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