<i>Instruções de Uso/Instructions for use/ Instrucciones de uso</i> RESINA ORTODÔNTICA ADITEK (LIGHTBOND) ADITEK ORTHODONTIC RESIN (LIGHTBOND) RESINA DE ORTODONCIA ADITEK (LIGHTBOND)	CÓDIGO/CODE/ CÓDIGO	IFU24
	SETOR/SECTOR/ SECTOR	Desenvolvimento/ Development/ Desarrollo
	PÁGINA/PAGE/ PÁGINA	1 de 15
	REVISÃO/REVIEW/ REVISIÓN	01

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Not applicable.

4. Usage Indication:

Aditek Orthodontic Resin is indicated for bonding brackets, tubes and accessories to teeth during orthodontic treatment.

The purpose of Aditek Orthodontic Resin is to promote the chemical and mechanical adhesion of the pieces to the respective teeth, which will serve as a support point for the application of other products used during orthodontic treatment.

NOTE: The choice of material to be used is made by the qualified professional, based on their knowledge acquired during their training together with the diagnosis and treatment plan for each patient.

5. Users:

The application of Aditek Orthodontic Resin is restricted to dental professionals. The use by people outside this context may cause harm to the patient.

6. Performance Features:

Aditek Orthodontic Resin is based on the fixation of brackets , tubes and accessories, promoting chemical and mechanical adhesion of the pieces to the tooth.

<p><i>Instruções de Uso/Instructions for use/ Instrucciones de uso</i></p> <p>RESINA ORTODÔNTICA ADITEK (LIGHTBOND) ADITEK ORTHODONTIC RESIN (LIGHTBOND) RESINA DE ORTODONCIA ADITEK (LIGHTBOND)</p>	CÓDIGO/CODE/ CÓDIGO	IFU24
	SETOR/SECTOR/ SECTOR	Desenvolvimento/ Development/ Desarrollo
	PÁGINA/PAGE/ PÁGINA	7 de 15
	REVISÃO/REVIEW/ REVISIÓN	01

As the resin is used to fix the pieces to the tooth, all the action done on the bracket, tube or accessories is transferred directly to the tooth.

With this transfer of force and torque, the position of the teeth is changed.

7. Clinical Benefits:

The main benefit that the orthodontic treatment performed brings to the patient's life is the aesthetic improvement with the alignment and leveling of the teeth, leaving the teeth of the dental arches in the correct and ideal position for a good occlusion.

8. Warnings:

 Aditek products are designed and manufactured for single use only. Its reuse or reprocessing can cause cross-infection, decrease in bonding efficiency and loss of properties.

 If the packaging has been violated, do not use the product, avoiding risks related to the transmission of microorganisms or the possible falsification of products.

Store the product in a cool, dry place away from sunlight, keeping it in its own packaging until use.

The product must be applied by a dentist who must follow the instructions for use and guidelines specified by the manufacturer.

The composition of the product does not require specific procedures for its transport and does not present risks to people.

The professional must use good clinical practices for the disposal of medical products according to current health regulations.

Aditek do Brasil Ltda. is not responsible for the results obtained by applying the product without proper precautions or failure to observe the warnings and instructions for use.

In case of serious incidents with the product, Aditek do Brasil S.A. must be notified to the competent authority in the national scope in which the dentist is established.

9. Contraindications:

- Patients unable to cooperate with treatment,
- Illness or limitation that may make treatment difficult,
- Decalcification of existing tooth enamel,
- Existing periodontal complications.

Note: May cause allergic reactions in patients with a history of sensitivity to acrylic and methacrylic resins.

<i>Instruções de Uso/Instructions for use/ Instrucciones de uso</i> RESINA ORTODÔNTICA ADITEK (LIGHTBOND) ADITEK ORTHODONTIC RESIN (LIGHTBOND) RESINA DE ORTODONCIA ADITEK (LIGHTBOND)	CÓDIGO/CODE/ CÓDIGO	IFU24
	SETOR/SECTOR/ SECTOR	Desenvolvimento/ Development/ Desarrollo
	PÁGINA/PAGE/ PÁGINA	8 de 15
	REVISÃO/REVIEW/ REVISIÓN	01

In case of allergic reaction of the patient, discontinue use and replace with similar material.
If persistent reactions occur, even after discontinuation of use, medical advice should be sought.

NOTE: Aditek do Brasil S.A. supplies products only to dental professionals. The dentist is entirely responsible for identifying and/or possible contraindications that may prevent the use of the product.

10. Residual Risks:

Not applicable.

11. Unwanted Side Effects:

May cause allergic sensitivity in individuals predisposed to acrylates and methacrylates due to prolonged or continuous use.

Signs and Symptoms of Exposure: Skin and eye irritation.

12. Composition:

LightBond - Folder	
Ingredients	Chemistry Composition Limits (%)
Bis Glycidyl Dimethacrylate Resin A (BIS-GMA) Polymer matrix forming monomer	15%
Triethylene Glycol Dimethacrylate (TEGDMA) Resin Polymer Matrix Forming Monomer	7%
Particulate barium borosilicate glass, surface treated with methacryloxypropyltrimethoxy silane - Reinforcement or filler	77%
Camphorquinone / Photoinitiator	0.3%
Ethyl 4-Dimethylaminobenzene (DMPT) / Co-initiator	0.5%
Titanium Dioxide - Pigment	0.1%
Iron oxides - Pigment	0.1%

Composition of the bonding agent - Sealant	
Ingredients	Chemistry Composition Limits (%)
Hydroxyethyl Methacrylate	89%
Deionized water	10.2%
Camphorquinone – Photoinitiator	0.3%
Ethyl 4-Dimethylaminobenzene (DMPT) / Co-initiator	0.5%

Conditioning Agent - Phosphoric Acid - Gel	
Ingredients	Chemistry Composition Limits (%)
Phosphoric acid	37%
Water	63%

13. Validity:

1 year when stored in its own packaging.

<i>Instruções de Uso/Instructions for use/ Instrucciones de uso</i> RESINA ORTODÔNTICA ADITEK (LIGHTBOND) ADITEK ORTHODONTIC RESIN (LIGHTBOND) RESINA DE ORTODONCIA ADITEK (LIGHTBOND)	CÓDIGO/CODE/ CÓDIGO	IFU24
	SETOR/SECTOR/ SECTOR	Desenvolvimento/ Development/ Desarrollo
	PÁGINA/PAGE/ PÁGINA	9 de 15
	REVISÃO/REVIEW/ REVISIÓN	01

14. Symbols and Meanings:

	Manufacturer		Do not reuse the product
	Consult instructions for use		Lot number
	Attention		Product reference code
	Do not use if the package is damaged or violated		Expiration date
	Protect from light		Temperature limit

15. Patient Information:

The product must be prevented from coming into contact with the skin, eyes and mucous membranes.

The professional must correctly follow the instructions for using the product. This product should only be used on patients by a dentist.

It is important to be careful when chewing certain types of food, especially hard or crunchy foods, as they can cause the resin to detach. These foods must be cut before eating. Hard foods like popcorn and ice should be avoided.

Do not put objects in your mouth, such as pens and pencils, and be careful not to hit or bite the toothbrush structure when brushing your teeth, which can break or loosen the bracket when subjected to a prying movement.

In case of ingestion or aspiration, seek medical attention.

16. Usage Mode:

For the correct use of the product and to obtain satisfactory results in the orthodontic treatment, carefully follow the steps described below:

- Correctly prepare the tooth that will receive the bracket, tube or other orthodontic accessory, using oil-free and fluoride-free paste, taking care not to injure the gingival tissue, to avoid contamination with blood from the area to be bonded;
- Wash well and isolate the tooth with cotton, leaving the surface of the tooth completely clean and dry;
- Apply phosphoric acid on the surface of the tooth with a brush or disposable sponge, to promote the superficial decalcification of the same, leaving it to act for approximately 30 to 60 seconds. This procedure will increase the adhesion of the resin to the tooth surface. If the tooth is not white, reapply the conditioner for another 20 seconds;

<p><i>Instruções de Uso/Instructions for use/ Instrucciones de uso</i></p> <p>RESINA ORTODÔNTICA ADITEK (LIGHTBOND) ADITEK ORTHODONTIC RESIN (LIGHTBOND) RESINA DE ORTODONCIA ADITEK (LIGHTBOND)</p>	CÓDIGO/CODE/ CÓDIGO	IFU24
	SETOR/SECTOR/ SECTOR	Desenvolvimento/ Development/ Desarrollo
	PÁGINA/PAGE/ PÁGINA	10 de 15
	REVISÃO/REVIEW/ REVISIÓN	01

- Then wash the tooth eliminating all acid residue;
- Isolate the tooth again with cotton, leaving the surface of the tooth completely clean and dry;
- Apply a thin layer of the liquid (bonding agent) to the tooth, using a disposable brush or sponge;
- Then, carefully apply the resin paste to the base of the piece, using a plastic or disposable wooden spatula;
- Position the piece correctly on the tooth with light pressure. Remove excess resin around the perimeter of the base of the piece;
- Apply the curing light on two sides (mesial - distal or incisal - gingival) for 20 to 40 seconds, using specific curing light equipment with a wavelength between 400 and 500 nanometers (blue light) with cycles of approximately 20 seconds.

NOTE: The use of protective eyewear is required by the professional and the patient. When using ceramic or plastic brackets, photopolymerization can be done directly on the bracket . Plastic brackets may need special primer on the base before bonding (not included in the kit).

- Remove the insulation and wash with plenty of water;
- Wait at least 5 minutes after gluing the last piece to place the arches and accessories;
- After being properly fixed, the pieces may receive the application of other products and accessories that complement the device in the orthodontic treatment, such as wires, elastic bands and arches, for the correction of dysfunctions and malocclusions of the dental arches.

16.1. Useful Tips:

- The intensity of the polymerization light may vary. Less intense light may require a longer curing time.
- The dispensed primer and paste should not be exposed directly to bright light for an extended period of time.
- On plastic brackets, use only primer specifically recommended for this purpose.
- When you prefer a firmer consistency of the resin, store it in the fridge and take it out only when using.

16.2. Bracket removal :

Brackets must be removed with suitable pliers, using minimal force. Avoid shearing or twisting forces. Greater care should be taken when removing ceramic brackets, as adhesion to tooth enamel is greater. Usually most of the adhesive is removed with the bracket. The remaining adhesive must be removed with appropriate instruments, or with a finishing bur.

<i>Instruções de Uso/Instructions for use/ Instrucciones de uso</i> RESINA ORTODÔNTICA ADITEK (LIGHTBOND) ADITEK ORTHODONTIC RESIN (LIGHTBOND) RESINA DE ORTODONCIA ADITEK (LIGHTBOND)	CÓDIGO/CODE/ CÓDIGO	IFU24
	SETOR/SECTOR/ SECTOR	Desenvolvimento/ Development/ Desarrollo
	PÁGINA/PAGE/ PÁGINA	11 de 15
	REVISÃO/REVIEW/ REVISIÓN	01

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