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1. Product's name:

1.1. Technical Name: Plate for Making Aligner Trays

1.2. Comercial Name: Memoflex Copolyester

1.3. Business Model: Memoflex Copolyester Plate 0.30mm, 0.40mm, 0.50mm, 0.60mm, 0.75mm, 0.76mm, 0.80mm, 1.00mm.

1.4. Anvisa Notification: 10331430037

2. Manufacturer:

Aditek do Brasil S.A.

Address: Cesário Mota, 14.

CEP: 14140-000 - Cravinhos - SP - Brazil

Phone: +55 (16) 3951-9355 / 0800-0159520

CNPJ 64.602.097/0001-95

Technical Manager: José Arantes Nogueira Neto - CREA: 601037885

www.aditek.com.br

3. Certifying entity:

Not applicable.

4. Indication of Use:

Memoflex Copolyester plates are indicated for making dental trays and, according to their characteristics and thickness, for different dental applications such as orthodontic aligners, occlusion plate, tray for indirect bonding of brackets, dental whitening plate and orthodontic retaining plates.

5. Users:

The use of plates for making trays is restricted to professionals in the field of dentistry. The use by people outside this context may cause harm to the patient. The user undertakes to verify his responsibility for the suitability of the product for the intended use, especially in the case of fields of use not expressly indicated in the instructions for use.

6. Performance Characteristics:

- Resistant to shocks and drops, even in small thicknesses;
- Optimum maintenance of the shape obtained after vacuum thermoforming (dimensional stability);
- Inert to the action of most chemicals;
- Excellent transparency;
- It practically does not cause allergic reactions or irritations in the tissues in contact;
- Hygroscopic, that is, it easily absorbs moisture.

7. Clinical Benefits:

Thermoformed polymer: Fidelity in the reproduction of the arcades.

Excellent translucency: Aesthetics during use.

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8. Warnings:

 Aditek products are designed and manufactured for single use only. Its reuse or reprocessing can cause cross infection and decreased efficiency.

 If the packaging has been violated, do not use the product, avoiding risks related to the transmission of microorganisms or possible counterfeiting of products.

Store the product in a place at room temperature, dry and away from sunlight, keeping it in its own packaging, so as not to change its shape.

This product must be transported in its own packaging, must not be folded or crushed and protected from moisture.

The product is supplied in a non-sterile condition.

The professional must use good clinical practices for the disposal of medical products in accordance with current health regulations.

Aditek do Brasil S.A. is not responsible for the results obtained by applying the product without the due precautions or failure to observe the warnings and instructions for use.

In the event of serious incidents with the product, Aditek do Brasil and the competent national authority must be reported.

9. Contraindications:

Allergic reactions are unlikely but cannot be excluded and should not be used by people allergic to plastic.

10. Residual Risks:

Not applicable.

11. Unwanted Side Effects:

Not applicable.

12. Composition:

Copolyester
Chemical composition limits (%)
100%

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14. Symbols and Meanings:

	Manufacturer		Product reference code
	Consult instructions for use		Expiration date
	Warning		Do not reuse the product
	Do not use if package is damaged or tampered with		Lot number
	Keep away from sunlight		

15. Patient Information:

During the plastification process, do not neglect the plastifying equipment. The material in the plates heats up and softens quickly, which can become unusable and even cause damage to the plasticizing equipment if it softens too much.

When working with any model of plate for making trays, during making and finishing the tray, always use a mask for your safety.

Persons not used to vacuum thermoforming procedures must read the instructions for use provided by the manufacturer before using the equipment.

The responsibility for the clinical treatment, technique for the construction of the tray, side effects and the action expected from the use of trays generated by Memoflex Copolyester plates are the sole responsibility of the dentist and/or the professional responsible for their use.

The metallized and vacuum-sealed plastic packaging protects the plates from moisture and contaminants. The plates made of Copolyester are highly hygroscopic (high moisture absorption capacity). After opening the vacuum package, the plates that are not used immediately must remain stored in the same package so as not to absorb moisture and not to present numerous small bubbles during heating in the plasticizer.

16. Mode of use:

All trays for making trays are wrapped with insulating film and can be processed in a pressure or vacuum thermoforming unit.

If processed in a pressure thermoformer, this extremely thin film must be thermoformed together with the plate, on the side of the model, ensuring an improved transparency of the plate, without the need for additional insulation of the model surface. After thermoforming, the film must be removed.

If processed in a vacuum thermoformer, the insulation film must be removed prior to thermoforming.

The setting of the heating and cooling time depends on the thickness of the thermoforming sheet and can be programmed directly in the thermoforming unit.

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16.1. Ideal Plasticization Point:

Due to the existence of several types of plates with different thicknesses, the ideal plastification point is identified by alterations in the plate, thus occurring different heating times, which must be evaluated and calibrated by the professional. In general, the sweet spot is observed by the change in texture and surface color (gloss). For less translucent boards, the board will look shiny and transparent. For thicker or opaque boards the ideal point is to drain the board by 10 to 12 mm.

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