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1. Product's Name:

1.1. Technical Name: Instrument for Dentistry

1.2. Commercial name: Turbo Bite

1.3. Commercial Model: 01 - Turbo Bite (1 Cable and 4 Tips)

1.4. Anvisa Notification: 10331430039

2. Manufacturer:

Aditek do Brasil S.A.

Address: Cesário Mota, 14.

CEP: 14140-000 - Cravinhos - SP - Brazil

Phone: +55 (16) 3951-9355 / 0800-0159520

CNPJ 64.602.097/0001-95

Technical Responsible: José Arantes Nogueira Neto - CREA: 601037885

www.aditek.com.br

3. Certifying entity:

Not applicable.

4. Usage Indication:

Bite Turbo devices are indicated to assist in the manufacture of a resin artifact that will disarticulate the occlusion, unlocking the teeth so that movement occurs without occlusal interference.

In cases of deep bite, the artifact is indicated to open the bite and to prevent and/or eliminate bracket interference.

In cases of class III pseudos, it is indicated to uncross the anterior bite, applying the resin to the lingual of the lower incisors.

NOTE: The choice of light-curing resin to be used is made by the qualified professional, based on their knowledge acquired during their training together with the diagnosis and treatment plan for each patient.

5. Users:

The use of Turbo Bite devices is restricted to dental professionals. The use by people outside this context may cause harm to the patient.

6. Performance Features:

The flat part of the mini mold must be facing the incisal of the upper tooth, so that the lower incisors occlude over the flat part of the artifact made by Bite Turbo, when closing the mouth, disoccluding the teeth.

The surface tilt angle can be changed back and forth by reversing the mini mold position on the positioning handle or using a diamond bur to reshape the resin surface for slight jaw repositioning.

The same procedure is indicated for pseudo class III cases. Use the smallest size lingual of the lower incisors, with the flat part facing the incisal.

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7. Clinical Benefits:

The clinical benefits that these systems bring to the patient's life are the aesthetic improvement with the alignment and leveling of the teeth, more efficient masticatory function facilitating swallowing, better digestion and prevention of occlusal problems.

8. Warnings:

Turbo Bite should help the molding of the resin to be bonded, in most cases, on the lingual phase of the maxillary central incisors, over the cingulate region.

These artifacts made with Turbo Bite are temporary and should remain in the mouth for a maximum of 5 to 7 months.

Inform the patient that the teeth will not touch after bonding the resin.

Instruct the patient to eat soft foods at first, until he adapts to the new situation and returns to eating normally.

Avoid pizza, chewing gum and sticky foods that can dislodge the resin.

If you notice wear on the incisal edge of the opposing teeth, remove the resins immediately.

If the packaging has been violated, do not use the product, avoiding risks related to possible falsification of products.

Store the product in a cool, dry place away from sunlight, keeping it in its own packaging until use.

The composition of the product does not require specific procedures for its transport and does not present risks to people.

The product is supplied in a non-sterile condition.

The professional must use good clinical practices for the disposal of medical products according to current health regulations.

Aditek do Brasil S.A. is not responsible for the results obtained by the application of the product, being the responsibility of the qualified professional for the correct use.

In case of serious incidents with the product, Aditek do Brasil must be notified to the competent authority in the national scope where the dentist is established.

9. Contraindications:

The product has no contraindications when used for its intended purpose, following its instructions for use.

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It is the sole responsibility of the qualified professional to identify and/or possible contraindications that may prevent the use of the product.

10. Residual Risks:

Not applicable.

11. Unwanted Side Effects:

Not applicable.

12. Composition:

Turbo Bite - Cable
Chemical composition limits (%)
Polypropylene (PP)
100%

Turbo Bite - Tip
Chemical composition limits (%)
elastomer
100%

13. Symbols and Meanings:

	Manufacturer		Product reference code
	Consult instructions for use		Lot number
	Heads up		Do not use if package is damaged or tampered with

14. Patient Information:

Turbo Bite devices do not pose any danger to their users and do not have restrictions on use as long as their limitations of use and their original purpose are observed, that is, the use and application of the instrument in situations not indicated for the instrument are not the responsibility of the manufacturer but from misuse by a qualified professional.

The Product is presented in PVC, plastic or cardboard packaging with an EVA cradle. Composed of 4 interchangeable tips.

The Turbo Bite has a cable length of 126,760mm with a diameter of 7,100mm. The tip has a diameter of 3,750 mm and a height of 3,780 mm.

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15. Usage Mode:

- 1- Select the most suitable size of the Turbo Bite mold from the 4 size options. Make sure that the size of the mold chosen will not allow the patient to accommodate the mandible behind the Bite Turbo after bonding, nullifying the effect of disarticulation.
- 2- Clean the lingual surface of the teeth with pumice stone.
- 3- Make the acid attack, according to the manufacturer's instructions.
- 4- Isolate, dry and apply the sealant.
- 5- Position the mini silicone mold on the instrument handle. The curved part of the mini mold should be positioned to follow the lingual gingival contour of the teeth and the flat part parallel to the incisal edge of the teeth.
- 6- Isolate the interior of the mini silicone mold cavity with Vaseline. Do not put Vaseline on the edges of the mini mold, only inside the cavity.
- 7- Fill the mini mold with light-curing resin and place it over the cingulate of the Upper Central Incisor, using the straight part of the mini mold as a reference, aligning it parallel to the incisal edge of the tooth.
- 8- Put the light to light cure.
- 9- Remove the mini mold carefully, parallel to the flat surface of the Turbo Bite.
- 10- Make contact adjustments on the manufactured artifacts, using carbon for marking and a diamond cutter for wear.

After use, the Turbo Bite device must be washed and disinfected.
It must not be autoclaved.

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