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1. Product Name:

1.1. Technical Name: Odontologic Acrylic Resins

1.2. Commercial Name: Self BioPrint

1.3. Commercial Model: Self BioPrint Guide; Self BioPrint IBT; Self BioPrint Splint Soft; Self BioPrint Splint Hard

1.4. Anvisa Notification: 10331439004

2. Manufacturer:

Aditek do Brasil S.A.

Address: Cesário Mota, 14.

CEP: 14140-000 - Cravinhos – SP - Brazil

Phone: +55 (16) 3951-9355 / 0800-0159520

CNPJ 64.602.097/0001-95

Technical Manager: José Arantes Nogueira Neto - CREA: 601037885

www.aditek.com.br

3. Certifying entity:

Not applicable.

4. Indication of use:

The Self BioPrint Odontologic Acrylic Resins are photo-activated resins curable by UV light, exclusive for use in 3D printers, which are indicated for the following applications:

- Self BioPrint Guide is indicated for making surgical guides, orthognathic surgery guides and implant positioning guides.
- Self BioPrint IBT is indicated for making guides or trays for indirect bonding of brackets.
- Self BioPrint Soft and Self BioPrint Hard are indicated for making myo-relaxing plates or bruxism plates.

The commercial models of this family are constituted of the same base of chemical components with variations of their concentrations to meet the specificity of each application.

NOTE: The choice of Self BioPrint acrylic resin to be used is made by the qualified professional, based on their knowledge acquired during their training in conjunction with the diagnosis and treatment plan for each patient.

5. Users:

The use of Self BioPrint acrylic resin is restricted to qualified professionals. Use by persons outside this context may result in damage to the user of the final product.

6. Performance Characteristics:

The working principle of Self BioPrint Acrylic Odontologic Resins in their process of obtaining each final product is by photo polymerization of the oligomers and monomers directed by ultraviolet light in the 3D printer, the format of the file previously established according to the

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patient's data and the use of each final product. The mechanism of action is by hardening the oligomers and monomers obeying the principle of chemical photo polymerization under UV light.

7. Clinical Benefit:

Not applicable.

8. Warnings:

If the packaging has been violated, do not use the product, avoiding risks related to the transmission of microorganisms or the possible counterfeiting of products.

The product must be used by a qualified professional who must follow the instructions for use and guidelines specified by the manufacturer.

The product must be transported in its original packaging, room temperature, low humidity and not be exposed to the sun or UV light. Must be stored and maintained at room temperature and low humidity and not be exposed to sun or UV light.

Heavy objects should not be deposited on the product packaging.

The professional must use good clinical practices for the disposal of medical products according to current sanitary standards to avoid contamination of land, plumbing and surface water.

The product can be toxic to aquatic organisms and can cause long-term adverse side effects on the aquatic environment (H413).

Aditek do Brasil S.A. is not responsible for the results obtained by the application of the product without due precautions or failure to observe the warnings and instructions for use.

Keep out of reach of children and animals.

Product supplied in a non-sterile condition and must be sterilized before use according to the method indicated by the manufacturer.

In the event of serious incidents with the product, Aditek do Brasil shall be notified to the competent national authority where the dentist is established.

9. Contraindications:

The product does not have contraindications when used for its intended purpose, following its instructions for use.

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It is the sole responsibility of the qualified professional to identify and/ or possible contraindication that may prevent the use of the product.

10. Residual Risk:

Not applicable.

11. Undesirable Side-effects:

Not applicable.

12. Composition:

Model	Composition
Self BioPrint Guide	Acrylic Monomers
	Acrylic Oligomers
	Photoinitiator
	4-methoxyphenol
	Pigments
Self BioPrint IBT	Acrylic Monomers
	Acrylic Oligomers
	Photoinitiator
	4-methoxyphenol
	Pigments
Self BioPrint Splint Soft	Acrylic Monomers
	Acrylic Oligomers
	Photoinitiator
	4-methoxyphenol
	Pigments

Model	Composition
Self BioPrint Splint Hard	Acrylic Monomers
	Acrylic Oligomers
	Photoinitiator
	4-methoxyphenol
	Pigments

13. Shelf Life:

 1 year

NOTE: The product is valid for 3 months after opening the bottle, if stored under appropriate conditions, at room temperature and not exposed to the sun or UV light or heat sources.

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14. Symbols and Meanings

	Manufacturer		Product Reference Code
	Consult instructions for use		Lot Number
	Attention		Shelf Life
	Do not use if package is damaged or tampered with		Keep away from sunlight

15. Information to patient:

Self BioPrint Odontologic Acrylic Resins are presented in high density Polyethylene bottles of 250, 500 and 1000g with liquid content.

Each package is provided with a self-adhesive label containing the product information.

It is a translucent and viscous liquid, with density between 300 and 1200 Cps, density between 0.9 and 1.2 g/cm³ and has a slight characteristic odor.

The product offers no danger to its users and has no restrictions on use provided that observed its limitations of use and its original purpose, that is, the use and application of the product in situations not indicated for the same are not the responsibility of the manufacturer but the misuse by the qualified professional.

It should not be used if the expiry date is expired.

The product must be handled with gloves, goggles and preferably in an airy environment.

The product can be irritating to the skin, eyes and respiratory system.

- In case of skin contact:
Dry with absorbent paper or towel;
Wash immediately with plenty of soap and water;
Do not use organic solvents.
In the case of dermatitis, seek medical help.
- In case of eye contact:
Rinse thoroughly with running water for 20 to 30 minutes.
Refer to the doctor if necessary.

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- In case of inhalation:
Go to a cool and airy place if there is strong inhalation of the product.
Refer to the doctor if necessary.
- In case of ingestion:
Do not cause vomiting.
The affected person should drink from 0.5L to 0.8L of water and seek medical help.

16. Mode of use:

- 01- The ideal working temperature of the resin should be from 20 C to 31 C;
- 02- The product must be handled in an environment free of UV light, preferably in a room with warm yellow or white light;
- 03- Use gloves, goggles, and an apron or lab coat to manipulate the product;
- 04- Check that the printer is calibrated with the print profile indicated by the manufacturer for the respective resin, so that the printer can be adjusted and print correctly.
- 05- The product must be shaken before use, to ensure the homogenization of its components;
- 06- Must be placed carefully in the printer tank, respecting the volume indicated by the manufacturer of each printer;
- 07- After the printing is finished, the product must be removed carefully from the printer, as it has not yet gone through the post-cure and still this fragile;
- 08- It should be washed with isopropyl alcohol or ethyl under mechanical agitation or cavitation (use of agitation with ultrasound), preferably in two consecutive baths, to ensure that there is no liquid resin in the final piece, each bath can be performed from 3 to 5 minutes, important not to exceed this time not to weaken the piece;
- 09- The piece must be dried with the use of compressed air or by resting on a paper towel;
- 10- If it has been used in the process, print media in the file must be removed after drying, removal at this stage is easier and does not weaken the final part;
- 11- After complete drying, it must be placed in a UV light chamber, for the final post-curing of 30 to 90 min;
- 12- After finishing the process the liquid resin should be stored in its bottle, do not leave the resin in the printing tank.

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