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## 1. Product's Name:

**1.1. Technical Name:** Odontologic Varnishes

**1.2. Commercial Name:** Horus BioPrint Glaze

**1.3. Commercial Model:** 01-Horus BioPrint Glaze - 1 amber glass bottle with brush lid, containing 10 ml of odontologic varnish; 02-Horus BioPrint Glaze - 1 amber glass bottle with brush lid, containing 20 ml of odontologic varnish; 03-Horus BioPrint Glaze - 1 amber glass bottle with brush lid, containing 30 ml of odontologic varnish; 04-Horus BioPrint Glaze - 1 amber glass bottle with brush lid, containing 40 ml of odontologic varnish; 05-Horus BioPrint Glaze - 1 amber glass bottle with screw cap, containing 10 ml of odontologic varnish; 06-Horus BioPrint Glaze - 1 amber glass bottle with screw cap, containing 20 ml of odontologic varnish; 07-Horus BioPrint Glaze - 1 amber glass bottle with screw cap, containing 30 ml of odontologic varnish; 08-Horus BioPrint Glaze - 1 amber glass bottle with screw cap, containing 40 ml of odontologic varnish; 09-Horus BioPrint Glaze - 1 amber glass bottle with dropper cap, containing 10 ml of odontologic varnish; 10-Horus BioPrint Glaze - 1 amber glass bottle with dropper cap, containing 20 ml of odontologic varnish; 11-Horus BioPrint Glaze - 1 amber glass bottle with dropper cap, containing 30 ml of odontologic varnish; 12-Horus BioPrint Glaze - 1 amber glass bottle with dropper cap containing 40 ml of odontologic varnish..

**1.4. Anvisa Notification:** 10331439006

## 2. Manufacturer:

Aditek do Brasil S.A.

Adress: Rua Cesário Motta, 14

CEP: 14140-000 - Cravinhos – SP - Brasil

Phone: +55 (16) 3951-9355 / 0800-0159520

CNPJ 64.602.097/0001-95

Technical Manager: José Vicente Junior - CREA: 5069373070

Site: www.aditek.com.br

## 3. Certifying Entity:

Not applicable.

## 4. Indication of Use

The product Horus BioPrint Glaze Odontological Varnish is a varnish used as a glaze (varnish) in acrylic resins of total, partial and provisional dentures, indicated to seal eventual porosities and to give the desired brightness to the material. This seal stabilizes the surface, is abrasion resistant and provides a high gloss.

## 5. Users:

The varnish is only used in a odontologic office by a trained professional.

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## 6. Performance Features:

The working principle of Horus BioPrint Glaze Varnish is based on the photopolymerization of oligomers and monomers directed by ultraviolet light, after its application by brush on the surface of the odontologic part obtained in a 3D printer. The mechanism of action is due to the hardening of oligomers and monomers following the principle of chemical photopolymerization under ultraviolet light.

## 7. Clinical Benefits:

The Horus BioPrint Glaze Odontologic Varnish has the clinical benefit of sealing porosities in total, partial and provisional acrylic resin dentures, providing high gloss on the parts.

## 8. Warnings:

If the packaging has been violated, do not use the product, avoiding risks related to the transmission of microorganisms or the possible falsification of products.

The product must be used by a qualified professional who must follow the instructions for use and guidelines specified by the manufacturer.

The product must be transported in its original packaging, at room temperature, low humidity and not exposed to sunlight or UV light. It should be stored and kept at room temperature and low humidity and not be exposed to sunlight or UV light.

Heavy objects must not be placed on the product packaging.

The professional must use good clinical practices for the disposal of medical products in accordance with current health standards to avoid contamination of the earth, plumbing and surface water. Aditek do Brasil S.A. is not responsible for the results obtained by applying the product without proper precautions or failure to observe the warnings and instructions for use.

The product must not be used if the expiry date has expired;

The product must be manipulated with gloves, goggles and preferably in a ventilated environment;

Keep out of reach of children and animals.

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#### 10. Residual Risks:

Avoid contamination of land, pipes and surface water;  
 Toxic to aquatic organisms, may cause long-term adverse side effects in the aquatic environment.

#### 11. Undesirable Side Effects:

Not applicable.

#### 12. Composition:

Chemical Name	Percentage by Weight (%)
Acrylated Monomers	> 96%
Pigments	≤ 1%
Photoinitiators	< 3%

#### 13. Validity:

2 years after manufacture if the bottle is unopened and 3 months after opening the bottle, protected from light, humidity and in a ventilated place.

#### 14. Symbols and Meanings

	Manufacturer		Do not reuse the product
	Consult instructions for use		Lot number
	Attention		Expiration date
	Do not use if package is damaged or tampered with.		Product reference code
	Protect from light		

#### 15. Patient Information:

The professional must correctly follow the instructions for use of the product. This product should only be used on patients by a dentist.

The product can be irritating to the skin, eyes and respiratory system;

In case of contact with the skin: dry with absorbent paper or a towel. Wash immediately with plenty of soap and water. Do not use organic solvents. In case of dermatitis, seek medical assistance;

In case of contact with eyes: rinse with plenty of running water for 20 to 30 minutes, refer to a doctor if necessary;

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In case of inhalation: move to a cool and airy place if there is strong inhalation of the product, refer to a doctor if necessary;

**16. How to Use:**

- 1 - Apply a thin uniform layer using a soft brush in one direction avoiding bubbles.
- 2 - Leave to act for 15 seconds before curing.
- 3 - Polymerize in UVA light apparatus for 10 to 15 minutes.

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