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1. Product Name:

1.1. Technical Name: Orthodontic Bracket.

1.2. Trade Name: Aditek Self-Ligating Orthodontic Tube.

1.3. Commercial Model: EasyClip+ Passive Weld Self-Ligating Single Tube; EasyClip+ Passive Glue Self-Ligating Single Tube; EasyClip+ Interactive Welding Self-Ligating Single Tube; EasyClip+ Interactive Glue Self-Ligating Single Tube.

1.4. Anvisa Notification: 10331430016

2. Manufacturer:

Aditek do Brasil S.A.

Address: Cesário Mota, 14.

Zip Code: 14140-000 - Cravinhos – SP - Brazil

Phone: +55 (16) 3951-9355 / 0800-0159520

CNPJ 64.602.097/0001-95

Technical Responsible: José Vicente Junior - CREA: 5069373070

Website: www.aditek.com.br

3. Certifying Body:

SIQ Ljubljana

SLO-1000 Ljubljana, Trzaska cesta 2

Phone: +386 1 4778 100

info@siq.si

www.siq.si

4. Indication for Use:

Self-Ligating Orthodontic Tubes are indicated as auxiliary components in orthodontic treatment. Its purpose is to serve as a fitting for pre-contoured wires and archwires and as a support point for orthodontic elastics or other auxiliary products used in orthodontic treatment.

The tubes are glued to the molar teeth or can be welded into the bands that will be fixed to the teeth, the force applied to them is transferred to the teeth, thus causing the change of their position, and consequently their correction to the desired position.

Self-ligating tubes have a clip that keeps the wire inside the slot, which can be passive (the clip is straight so that the bow is not constantly pressed into the slot) or interactive (the clip is curved, causing the bow to be pressed into the slot).

These products were developed for exclusive use inside the oral cavity, to assist in orthodontic treatment, and for human use.

They are indicated to treat patients who have malocclusion problems such as prognathism, retrognathism, deep bite, open bite, crowding, anterior crossbite, posterior crossbite, and diastemas.

NOTE: The choice of the material to be used is made by the qualified professional, based on his/her knowledge acquired during his/her training together with the diagnosis and treatment plan for each patient.

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5. Users:

The application is restricted to dental professionals. Use by people outside this context may cause harm to the patient.

6. Performance Characteristics:

Usually, an orthodontic treatment without interruption takes from 1 to 3 years, varying in each case according to the facial pattern, type of malocclusion, and oral and systemic conditions of the patient, and may extend according to the complexity of the case.

7. Clinical Benefits:

The main benefits that these systems bring to the patient's life are aesthetic improvement with the alignment and leveling of the teeth, better hygiene, more efficient chewing function facilitating swallowing, better digestion, prevention of occlusal problems, better speech, and breathing.

8. Warnings:

Product containing Chromium and Nickel.

Aditek's products are designed and manufactured for single use.

If the packaging has been tampered with, do not use the product.

Store the product in a cool, dry place away from sunlight, keeping it in its own packaging until use.

The product is not recommended in the MRI setting.

The product is safe for use in pregnant women.

The composition of the product does not require specific procedures for its transport and does not present risks to people.

The product is supplied in the non-sterile condition.

The professional must use good clinical practices for the disposal of medical products in accordance with current health regulations.

Aditek do Brasil S.A. is not responsible for the results obtained by the application of the product without due precautions or non-observance of warnings and instructions for use.

In the event of serious incidents with the product, Aditek do Brasil S.A. must be reported to the competent authority in the national scope in which the dentist is established.

9. Contraindications:

- Patients with poor oral hygiene.
- Patients with inability to cooperate with treatment.

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- Patients who have an allergy to chromium and/or nickel.
- Illness or limitation that may hinder treatment.
- Existing bone or root resorption.
- Descaling of existing tooth enamel.
- Existing periodontal complications.

NOTE: Aditek do Brasil S.A. supplies products only to dental professionals. It is the sole responsibility of the dentist to identify and/or possible contraindications that may prevent the use of the product.

10. Precautions:

The product may cause allergic reactions in patients susceptible to Chromium and Nickel allergy. In case of reaction, remove the device and advise the patient to seek medical attention.

The reuse or reprocessing of the product can cause cross-infection, decreased bonding efficiency, and loss of mechanical properties due to natural wear and forces applied in the removal of the product.

Avoid risks related to the transmission of microorganisms or possible counterfeiting by using only products that have their packaging sealed.

The product should be applied by a dentist who should follow the instructions for use and guidelines specified by the manufacturer.

In the MRI environment, the product will reflect the magnetic waves generated by the CT scanner, preventing, or masking the correct reading of the scan, so they should always be removed before MRIs of the head and neck.

Due to hormonal changes during pregnancy, a woman is more susceptible to periodontal problems. During orthodontic treatment, the difficulty of cleaning can aggravate these problems. It is important that the patient is referred for professional cleaning at regular intervals.

The procedure of bonding and removing the tubes on the teeth must be performed according to the manufacturer's instructions for use of the orthodontic adhesive purchased.

The product is supplied clean in a non-sterile condition, and it is at the discretion of the professional to perform the sterilization.

11. Residual Risks:

All finished orthodontic treatments are subject to recurrence. After treatment, control should be determined by the professional to avoid or minimize recurrence.

12. Undesirable Side Effects:

The use of the medical device can lead to discoloration or decalcification of teeth, bone and root resorption, periodontal complications, allergic reactions, difficulties in maintaining oral hygiene, damage to the oral mucosa and tongue, difficulty speaking or chewing, discomfort and pain.

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13. Composition:

AISI 17-4PH Stainless Steel – Tube Body								
Chemical Composition Limits (%)								
C	Mn	P	S	Oneself	Cr	Nu	Ass	Fe
≤ 0.070	≤ 1.00	≤ 0.040	≤ 0.030	≤ 1.00	15.0 to 17.5	3.0 to 5.0	3.0 to 5.0	Remainder
ASTM: A564								

Nickel-Titanium - Clip									
Chemical Composition Limits (%)									
Nu	C	Co	Ass	H	Fe	Nb	N+O	Cr	Thee
54.5 to 57.0	≤ 0.050	≤ 0.050	≤ 0.010	≤ 0.005	≤ 0.050	≤ 0.025	≤ 0.050	≤ 0.010	Remainder
ASTM: F2063									

14. Validity:

10 years.

15. Symbols and Meanings:

	Manufacturer		Authorized representative in the European Community
	See instructions for use		Complies with the European directive. XXXX represents the number of the certifying body
	Attention		Product contains chromium and nickel
	Do not use if package is damaged or tampered with		Do not reuse the product
	Product Reference Code		Batch number
	Expiration date		

16. Patient Information:

Product containing Chromium and Nickel in its composition, which can cause allergic reactions in people susceptible to these components.

The patient must perform good daily hygiene to avoid the installation of bacterial plaque and consequent periodontal diseases.

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It is important that the pregnant patient is referred for professional cleaning at regular intervals, due to the increased susceptibility to periodontal problems during pregnancy.

It's important to be careful when chewing certain types of foods, especially harder or crunchy ones, as they can peel off or break the tubes. These foods should be cut before eating. Hard foods, such as popcorn and ice, should be avoided.

Do not put objects in your mouth, such as pens and pencils, and be careful not to hit or bite the brush frame when brushing your teeth, which can break or loosen the tube when subjected to a lever.

In the event of a breakdown, seek a dentist.

In case of ingestion or aspiration, seek medical attention.

When participating in sports, always wear mouth guards, as impacts can damage the braces and injure the oral cavity.

The patient should follow the professional's instructions regarding restraint to avoid recurrence.

The product is not recommended in the MRI setting.

17. Mode of Use:

The steps below should only be performed by dental professionals.

The tubes should be glued directly to the surface of the teeth at their base, with the use of orthodontic adhesives, following the instructions for use by the adhesive manufacturer.

To open the clip: Place the instrument (explorer) in the clip hole, located below the slot, between the lower fins and perpendicular to the buccal surface of the tooth. Gently rotate the instrument to the right or left; if the workpiece has a hook, the rotation should be in the opposite direction to the hook. The clip can also be opened with a fine-tipped explorer with a gentle downward lever movement, placing the instrument inside the hole. The lower tube clips open to the gingival and the upper ones to the occlusal.

Position the selected arch according to the treatment to be performed.

Clip Closure: Using tweezers, place one end over the fins of the tube and the other end into the clip. Press the tweezers vertically. The first "click" unlocks the opening that was made in the clip. The second, stronger "click" will lock the clip, closing the slot.

After treatment, remove the tubes and orthodontic adhesive from the teeth using an appropriate instrument, according to the manufacturer's instructions for the use of the orthodontic adhesive purchased.

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